

# **Cell Counting using the Trypan Blue Exclusion Method**

# Learning outcomes



- Understand the basic structure of hemocytometer.
- Learn to load sample on hemocytometer.
- Defining Neubauer factor?
- Understand formulae for calculating cell viability and cell number in a sample.
- Calculate cell numbers using hemocytometer.
- Calculate cell viability in a sample after staining the cells with trypan blue.
- Practising focusing an inverted microscope to visualize mammalian cells.
- Identify the morphology of cells (sample given during practical)
- Estimating the confluency of cells in three different samples.

# Trypan Blue

- Trypan Blue Solution, 0.4%, is routinely used as a **cell stain** to assess cell viability using the **dye exclusion test**.
- This test is often performed while counting cells with the **hemocytometer** during routine **sub-culturing**, but can be performed any time **cell viability** needs to be determined quickly and accurately.
- The dye exclusion test is based upon the concept that **viable cells do not take up impermeable dyes** (like Trypan Blue), but dead cells are permeable and take up the dye.


# Hemocytometer

- The **hemocytometer** or haemocytometer is a device originally designed for the counting of blood cells. It is now also used to count other types of cells as well as other microscopic particles.
- The hemocytometer was invented by Louis-Charles Malassez and consists of a **thick glass microscope slide** with a **rectangular indentation** that creates a chamber. This chamber is engraved with a laser-etched grid of perpendicular lines. The device is carefully crafted so that the area bounded by the lines is known, and the **depth of the chamber is also known**. It is therefore possible to count the number of cells or particles in a specific volume of fluid, and thereby calculate the concentration of cells in the fluid overall.



# Neubauer – Hemocytometer

- In an improved Neubauer hemocytometer, the total number of cells per ml can be discovered by simply multiplying the total number of cells found in the hemocytometer grid (16 Squares) by  $10^4$ .




<u>Hemocytometer grid:</u>	
Red square	= $1 \text{ mm}^2$
Green square	= $0.0625 \text{ mm}^2$
Yellow square	= $0.04 \text{ mm}^2$
Blue square	= $0.0025 \text{ mm}^2$

*Count cells in 4 corner squares, on top and left border, take average*

# Neubauer Factor ( $10^4$ )

- With the coverslip on the chamber, the coverslip sits 0.1mm over the chamber. Therefore, the volume of each square, contained under the coverslip is:
- $1\text{mm} \times 1\text{mm} \times 0.1\text{mm} = 0.1\text{mm}^3 / 0.1\mu\text{l} / 10^{-4}\text{ ml}$   
**( $10\text{mm} = 1\text{cm}$ ;  $1\text{mm} = 0.1\text{cm}$ ;  $1\text{mm}^3 = 10^{-3}\text{ ml}$ ;  
 $0.1\text{mm}^3 = 10^{-4}\text{ ml}$  (the vol. of 1 square))**



# Neubauer Factor ( $10^4$ )

- If  $0.1 \mu\text{l}$  has = 1 cell
- Then  $1000 \mu\text{l}$  would have =  $\frac{1 \times 1000}{0.1} = 10000 = \mathbf{10^4}$  cells


So, for calculating cell number per mL, multiply the counted number of cells with **Neubauer Factor ( $10^4$ )**

# Cell Count Determination - Hemocytometer



Neubauer Hemocytometer

Cells in 16 squares counted



Light Microscope: Magnification 100 x

You got your cell suspension (10 ml) and want to know how many cells you might have in your cell suspension.

You add 1 Vol of cell suspension to 1 Vol of Trypan Blue solution.

**Trypan Blue** stains dead cells only.

You count cells in **16 squares**: e.g. **80 white** (unstained/white cells) And **20 blue cells** (stained cells). Makes a total of **100 cells**.

**Total Cell count determination (total cell concentration):**

$100 \times 2$  (due to 1:1 dilution) = 200

$200 \times 10^4$  (**Neubauer Factor**) =  $2 \times 10^6/\text{ml}$  (total cell concentration)

**Total Cell Number:**

**10 ml of cell solution:**  $10 \times 2 \times 10^6 = 2 \times 10^7$  cells  
(total cell number)


**Cell Viability (%):**  $\frac{\text{number of white cells} \times 100}{\text{total cell number}}$

$80 \times 100 / 100 = 80\%$  Cell Viability




# Inverted Microscope

- An **inverted microscope** is a microscope with its light source and condenser on top, whereas the objectives are below the stage pointing up. It was invented in 1850 by J. Lawrence Smith.
- Inverted microscopes are useful for observing living cells (or organisms) at the bottom of a tissue culture flask under **more natural conditions than on a glass slide**, as is the case with a conventional light microscope.



# Light Microscope



# Inverted Microscope



# Laboratory Practical

The purpose of this laboratory exercise is to acquaint you with some of the fundamentals of culturing cells *in vitro*.

- (1) What do cell cultures look like?
- (2) Where do the cells grow?
- (3) What do the cells in the culture look like?

In this exercise you will:

- (1) Observe the cell culture provided;
- (2) Count cells to determine how many cells were in the original plate;
- (3) Determine the viability of cells.

# MATERIALS

- Hemocytometer plus a supply of cover slips.
- 0.4% Trypan Blue stain (fresh & filtered) in phosphate buffered saline.
- Cell Suspension.
- Micropipettes and tips.
- Inverted microscope.

Lab Coats  
Gloves

## Workflow to be followed

**Visualize cells under microscope (record confluency, morphology, cell line name)**



**Prepare the hemocytometer**



**Prepare a 1:1 dilution of the cell suspension in trypan blue**



**Load the hemocytometer**



**Determine the number of cells (total, live and dead)**





**Calculations**



**Cleaning the hemocytometer**

# Cell Counting

- For purposes of this practical, you will only count cells touching the top and left border (count 'black' and not 'red' as shown in figure). You will **count the cells in four large corner squares** (i.e. those divided further into sixteen small squares).



# Calculations

## *Calculate the average #cells from chamber*

Average #

Live Cells 150

Average #

Dead Cells 20

Average #

Total Cells 170

## *Calculate the cell density or Total cell concentration in the original suspension (in cells/ml) is:*

# Cells/ml = (Average # total cells) X dilution factor X  $10^4$

Dilution Factor = Total Volume (Volume of sample + Volume of diluting liquid) / Volume of sample. The dilution factor for this example is 2 because 20  $\mu$ l of cell suspension was diluted with 20  $\mu$ l Trypan Blue.

# Cells/ml = 170 X 2 X  $10^4$

# Cells/ml =  $340 \times 10^4$  =  $3.4 \times 10^6$ /ml

## *Calculate cell viability*

Average # live cells / Average # total cells x 100

150 / 170 x 100

Cell viability = 88.23%

