

Experiment A3: Neutralisation Back Titration


Keelan Krinsky: 1634953

Full title: Analysis of Antacid tablets.

Demonstrator: Raymond

Group: 7

Date Performed: 04/10/2018



Aims

The main aim of this experiment was to establish the amount/mass of calcium (and magnesium) carbonate present in commercial antacid tablets.

Procedure

1. Gentle boiling of analyte solution to dissolve the antacid tablets was carried out for 10 minutes (5 minutes longer than the 5 minutes indicated in the laboratory manual)

Table 1: Volumes of 0.4997 standard HCl solution added in the neutralisation of Antacid analyte solution

	Tablet 1	Tablet 2	Tablet 3
Initial Reading/ml	37.74	37.59	36.40
Final Reading/ml	0.27	1.10	0.21
Volume Added/ml	37.47	36.49	36.19
hydrochloric acid added/mmol	18.72	18.23	18.08

Table 2: Volumes of 0.1049 standard NaOH solution added in the neutralisation of Antacid analyte solution, neutralised by excess strong acid

	Tablet 1	Tablet 2	Tablet 3
Initial Reading/ml	37.74	37.59	36.40
Final Reading/ml	0.27	1.10	0.21
Volume Added/ml	37.47	36.49	36.19
Amount of sodium hydroxide added/mmol	3.93	3.83	3.80
Amount of hydrochloric acid neutralised by tablet/mmol	14.79	14.40	14.28

Concentration of NaOH used = 0.1049M

Calcium carbonate content of antacid tablets.

CaCO_3 is assumed to be the only active weak base in the tablets provided

Average amount of hydrochloric acid neutralised per tablet = 14.49 mmol

Standard deviation in the measure of the amount of acid neutralised = 0.27 mmol

Relative standard error in the measure of the amount of acid neutralised = 1.84%

Mass of calcium carbonate per tablet = (Moles of calcium carbonate) \times (Molecular mass of calcium carbonate)
= 0.5 \times (Moles of HCl neutralised) (100.086 g \cdot mol⁻¹) (as HCl react 2:1 with CaCO_3)
= 0.5 \cdot 14.49 mmol \cdot 100.086 g \cdot mol⁻¹
= 725.1 mg

Questions


Why can the antacid tablets not be titrated directly?

The Antacid tablets contain only weak bases, in this case CaCO_3 the reaction between these weak bases and any strong acid which they were titrated against would not be fast enough for the practical determination


of an accurate endpoint in the titration, and hence accurate results overall. That is, the slow neutralisation would imply that the solution must be left to stand for long period after the addition of each drop of acid to determine if all of the $CaCO_3$ present has really been released.

Why does the carbon dioxide need to be driven off? Why would it cause problems with the titration?

The CO_2 must be driven off because, in an aqueous environment it also has acidic properties. Carbon dioxide reacts with water to form carbonic acid by the following equation:



The carbonic acid can then act to neutralise the strong base added in the back titration as follows:



This side reaction would lead to an inaccurately high amount of $NaOH$ added in neutralisation, and a subsequent underestimation of the Calcium carbonate content of the antacid tablets.

Tablet manufacturers are usually permitted to have an error of 7.5% in the amount of active ingredient in a tablet. Do all your three values for millimoles of HCl neutralized lie within $\pm 7.5\%$ of your mean value?

Tablet 1

$$\text{Percentage difference from mean} = |2.07\%| < 7.5\%$$

Tablet 2

$$\text{Percentage difference from mean} = |-0.62\%| < 7.5\%$$

Tablet 3

$$\text{Percentage difference from mean} = |-1.45\%| < 7.5\%$$

The percentage deviation for all of the tablets was well below the 7.5% indicated by the manufacturer.

Calculate the mass of $CaCO_3$ required to make up 500 mL of a solution containing 1000 ppm Ca.

$$\text{Moles of Ca in 500.0 mL of 1000 ppm Ca containing solution} = 1000mg \cdot L^{-1} \cdot 0.5000L = 500.0mg$$

$$\text{Moles of Ca used in preparation} = \frac{\text{Mass of Ca used}}{\text{Molar Mass of Ca}} = \frac{0.5000g}{40.078g \cdot mol^{-1}} = 0.012 mol$$

$$\text{Moles of } CaCO_3 \text{ required} = \text{Moles of Ca required} = 0.012 mol$$

$$\text{Mass of Moles of } CaCO_3 \text{ required} = (\text{Moles of } CaCO_3 \text{ required}) \times (\text{Molecular mass of } CaCO_3) = 0.012 mol \cdot 100.086g \cdot mol^{-1} = 1.249g$$