Experiment A2: Complexometric Titration

Keelan Krinsky: 1634953

Full title: Determination of Hardness of water

Demonstrator: Siya

Group: 7

Date Performed: 03/10/2018



Aims

- 1. Preparation of $\approx 0.005M$ EDTA standard
- 2. Determination of Calcium and Magnesium concentration in tap water by titration with EDTA, at pH 10 using Mercurochrome Black T indicator.
- 3. Determination of Calcium concentration in tap water by titration with EDTA standard in alkaline solution using demure indicator.
- 4. Determination of Magnesium concentration in tap water by difference calculation

Introduction

Tap water contains several mineral ions in reasonable high concentrations in particular mineral ions containing Calcium and magnesium, which enters the water bodies through the natural weathering of calcium and magnesium containing rocks (Geoffrey, 2007). Calcium and Magnesium concentrations form an important component in the monitoring of drinking water because they are closely related to the development of a number of serious diseases including heart disease and various forms of cancer (Leurs ,2010) Within water bodies calcium and magnesium typically exist in the form of suspended of dissolved salts such as carbonates. (McCreaon J. M. 1950). These salt forms are prone to precipitation, and as such as said to constitute the hardness of the water, although as Magnesium concentration is typically considerably lower than Calcium concentration hardness is often expressed merely as Calcium ion concentration. (Weingartner, 2006)

Calcium (and or Magnesium) ion concentration can be assessed by titration with a suitable chelating agent, such as for example Ethylenediaminetetraacetic Acid (EDTA). EDTA is a useful chelating agent as it will readily and rapidly react in a one to one ratio with metal cations to form reasonably stable complexes. However EDTA must be in its fully deprotonated form to form such complexes and as such a strongly basic environment must be maintained throughout the titration (Yoe, 1958).

The issue associated with EDTA however is that it is not specific to a given metal cation but will complex with all metal cations present in solution. Particularly in the case of determining Ca concentration any EDTA added will complex with Mg as well as Ca present in the analyte. This issue can be overcome by the use of an appropriate indicator. For example instead of an indicator such as Eriochrome Black which changes color only when EDTA is in excess and therefore available to complex with it, an indicator such as murexide forms a stable colored complex with Mg ions in the sample. Addition of EDTA first acts to chelate Ca, and then only acts to attack the stable magnesium complexes formed leading to a color change marking the endpoint of the titration. (Flaschka, 2013)

Results

Experimental Results

1) Preparation of EDTA standard solution.

Mass of EDTA used in standard= 9.9695g - 9.5254g = 0.4441g

Table 1: Volumes of EDTA standard solution titrated in the determination of Ca and Mg content of tap water

Titration Number	1	2	3
Final Reading/ml	14.29	27.52	13.30
Initial Reading/ml	0.25	14.30	0.18
Volume of EDTA solution delivered/ml	14.04	13.22	13.12

Table 2: Volumes of EDTA standard solution titrated in the determination of Ca content of tap water

Titration Number	1	2	3
Final Reading/ml	8.68	17.18	32.38
Initial Reading/ml	0.21	8.68	24.00
Volume of EDTA solution delivered/ml	8.47	8.50	8.38

- 2) Titration of Tap water with EDTA standard solution at pH 10 using Eriochrome Balck T indicator.
- 3) Titration of Tap water with EDTA standard solution under alkaline conditions using murexide indicator to determine Calcium content.

Calculations

Concentration of EDTA standard solution

Concentration of EDTA standard solution= $\frac{\frac{0.4441g}{372.24g \cdot mol^{-1}}}{0.2500L} = 0.0047722M$

Concentration of Mg and Ca in Tap water

Titration #1

Moles of EDTA titrated

Moles of EDTA titrated=0.067002mol

Moles of Ca and Mg in tap water aliquot

Moles of Ca and Mg =Moles of EDTA, as both Ca^{2+} and Mg^{2+} react one to one with EDTA when forming complex ions.

Moles of Ca and Mg in 100ml tap water aliquet = 0.067002mol

Concentration of Ca and Mg in tap water aliquut

Concentration of Ca and Mg in tap water aliqout= $\frac{0.067002mol}{0.1000L}$ =0.67002M

Titration #2

Moles of EDTA titrated

Moles of EDTA titrated=0.063088mol

Moles of Ca and Mg in tap water aliquot

Moles of Ca and Mg =Moles of EDTA, as both Ca^{2+} and Mg^{2+} react one to one with EDTA when forming complex ions.

Moles of Ca and Mg in 100ml tap water aliquet= 0.063088mol

Concentration of Ca and Mg in tap water aliquut

Concentration of Ca and Mg in tap water aliqout= $\frac{0.063088mol}{0.1000L}$ =0.63088M

Titration #3

Moles of EDTA titrated

Moles of EDTA titrated=0.062611mol

Moles of Ca and Mg in tap water aliquot

Moles of Ca and Mg =Moles of EDTA, as both Ca^{2+} and Mg^{2+} react one to one with EDTA when forming complex ions.

Moles of Ca and Mg in 100ml tap water aliquet = 0.062611mol

Average concentration of Mg and Ca in tap water.

Average concentration of Mg in tap water = 1.257M

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

Relative standard Error

Relative standard error for the concentration of Ca and Mg in tap water = 0.54%

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

Concentration of Ca and Mg in tap water aliquut

Concentration of Ca in tap water aliqout= $\frac{0.062611mol}{0.1000L}$ =0.62611M

Concentration of Ca in Tap water

Titration #1

Moles of EDTA titrated

Moles of EDTA titrated=0.04042mol

Moles of Ca in tap water aliquot

Moles of Ca =Moles of EDTA, as both Ca^{2+} and Mg^{2+} react one to one with EDTA when forming complex ions.

Moles of Ca in 100ml tap water aliquet= 0.04042mol

Concentration of Ca in tap water aliquut

Concentration of Ca in tap water aliqout= $\frac{0.04042mol}{0.1000L}$ =0.4042M

Titration #2

Moles of EDTA titrated

Moles of EDTA titrated=0.040564mol

Moles of Ca in tap water aliquot

Moles of Ca =Moles of EDTA, as both Ca^{2+} and Mg^{2+} react one to one with EDTA when forming complex ions.

Moles of Ca in 100ml tap water aliquet= 0.040564mol

Concentration of Ca in tap water aliquit

Concentration of Ca in tap water aliqout= $\frac{0.040564mol}{0.1000L}$ =0.40564M

Titration #3

Moles of EDTA titrated

Moles of EDTA titrated=0.039991mol

Moles of Ca in tap water aliquot

Moles of Ca =Moles of EDTA, as both Ca^{2+} and Mg^{2+} react one to one with EDTA when forming complex ions.

Moles of Ca in 100ml tap water aliquet= 0.039991mol

Concentration of Ca in tap water aliquit

Concentration of Ca in tap water aliqout= $\frac{0.039991mol}{0.1000L}$ =0.39991M

Average concentration of Ca in tap water.

Average concentration of Mg in tap water = 0.4028M

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

CaCO₃ Concentration in tap water

acMed<-c(cMed[2], cMed[3])

Concentration of $CaCO_3 = \frac{0.0406, 0.04mol \cdot 100.0869g \cdot mol^{-1}}{0.1000L} = 40.36mg \cdot L^{-1}$

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

Relative standard Error

Relative standard error for the concentration of Ca and Mg in tap water = 0.74%

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

Concentration of Mg in Tap water aliquot

Titration #1

Concentration of Mg= (Concentration of Mg and Ca)- (concentration of Ca)=0.67002 - 0.4042M = 0.26581M

Titration #2

Concentration of Mg = (Concentration of Mg and Ca)- (concentration of Ca)=0.63088-0.40564M=0.22525M

Titration #3

Concentration of Mg= (Concentration of Mg and Ca)- (concentration of Ca)=0.62611 - 0.39991M = 0.2262M

Average concentration of Mg in tap water.

Average concentration of Mg in tap water = 0.2257Mi

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

Relative standard Error

Relative standard error for the concentration of Ca and Mg in tap water = 0.3%

Titration 1 was excluded as a rough titration in this calculation, only the accurate titrations, titrations 2 and 3 where used.

Discussion

Although the true concentration of calcium within the sample is unknown and hence a definitive comment on the accuracy of the determination cannot be made, the concentration of Calcium (in the form of calcium carbonate) calculated fall within a reasonable range of expected values. The value of 40.35ppm calculated implies classification of the water a Soft (between 0 and 50 ppm $CaCO_3$), which is what would typically be expected of drinking water which has been treated for pallatability.

The relatively low relative standard error for the accurate titrations seen ($\leq 1\%$ in all cases) also serves to indicate that a reasonably high level of precision was achieved in the results. Although the low number of replicates implies less confidence that the high level of precision seen was not achieved by chance. The small degree of variation which was observed may be explainable by variations in the water samples used (different aliquots were drawn from different tabs), variations in the concentration of the EDTA standard due to insufficient mixing before titration, and difficulties in pinpointing the endpoint of titration, due to variations is the amount of indicator solution added.

In conclusion, the value obtain for calcium carbonate concentration in drinking water appears to be plausible and precise, identifying the sample in question as soft water.

References

McCreaon J. M. (1950), the Isotopic Chemistry of Carbonates and a Temperature Scale, *The Journal of Chemical Physics* 18:6, 849-857 Geoffrey M. G., (2007), Geomycology: biogeochemical transformations of rocks, minerals, metals and radionuclides by fungi, bioweathering and bioremediation, Mycological Research, 111.1: 3-49

Weingartner H (2006), "Water", Encyclopedia of Industrial Chemistry, Wiley-VCH, Weinheim

Leurs L. J., Schouten L. J., Mons M. N., Goldbohm R. A. Van den Brandt P. A. (2010), Relationship between Tap Water Hardness, Magnesium, and Calcium Concentration and Mortality due to Ischemic Heart Disease or Stroke in the Netherlands, *Environmental Health Perspectives*, 118, 3:414-420

Yoe, H. J., (1958), The Analytical Uses of Ethylenediaminetetraacetic Acid, J. Am. Chem. Soc, 80, 10:2600–2600

Flaschka, H. A. (2013). EDTA titrations: an introduction to theory and practice. Elsevier.